

Saphir 320/330 Series Grinding & Polishing Machines

10
inch



Hardy. Practical.

The **Saphir 320** (single wheel) and **Saphir 330** (dual wheel) manual preparation machines are grinder/polishers built for heavy usage and high-speed preparation. Available in either 200 mm (8") or 250 mm (10") wheel variable speed formats, the Saphir 320/330 series can be tailored for any application.

A solid, welded aluminum plate housing provides stability and corrosion resistance with a powder coated surface. All electrical components, spindles and moving parts are modular for easy cleaning and maintenance.

The Saphir 320/330 series is an excellent choice for coarse grinding.



Saphir 320/330 Variable Speed Grinding & Polishing Systems

320 Series - Single Wheel

Wheel / Disc Size:
200, 250 mm / 8, 10"

Speed:
50-600 rpm, Variable

Power:
110 V / 60 Hz, 1 Phase
220 V / 60 Hz, 1 Phase

Dimensions W x H x D:
16 x 11 x 25"

Weight: 49 lbs

330 Series - Dual Wheel

Wheel / Disc Size:
200, 250 mm / 8, 10"

Speed:
50-600 rpm, Variable

Power:
110 V / 60 Hz, 1 Phase
220 V / 60 Hz, 1 Phase

Dimensions W x H x D:
30 x 11 x 25"

Weight: 82 lbs



Wide open **Polishing and Grinding Wheel** with an extendable spigot for easy cleaning



Simple **Controls** and high torque motors will provide years of dependable service



Optional **Sedimentation Tank** for removing grinding waste

