

# Topas 130 Dosing System

Automatic



## Simple. Integrated.

The **Topas 130 Automatic Dosing System** makes dosing diamond abrasive easy. With space for three diamond suspensions and one lubricant, each solution can be adjusted to dose the perfect amount

for proper preparation of your samples.

While the **Topas 130** is designed to integrate flawlessly with ATM's line of grinder/polishers, it can also be used with any other brand.

## Topas 130 Automatic Dosing System

Variable Dosing:  
Manual or Automatic  
using Peristaltic Pumps

Dosing Intervals:  
2 seconds - 5 Minutes

Number of nozzles:  
3 Suspensions  
1 Lubricant

Power:  
110 V / 60 Hz, 1 Phase

Dimensions W x H x D:  
7 x 10 x 8"

Weight: 5 lbs



Individually control the  
precise, variable dosing of  
suspensions and lubricants  
with two second to five  
minute dosing intervals



Exact droplet dosing with  
reverse rinsing function  
keeps lines clear and ready  
for use



Mager's line of high  
performance diamond  
suspensions are made  
specifically for manual and  
automatic dosing systems

